



Hand in Hand for a Better Life!



**DEGANIA™  
LOW PROFILE  
GASTROSTOMY TUBES  
AND ACCESSORIES**

## INTRODUCTION TO ENTERAL FEEDING

Low Profile Gastrostomy Tubes ("G-buttons") are used by physicians in pediatric and adult patients with feeding/swallowing disorders. Typically, Low Profile Gastrostomy tubes are used when enteral feeding is expected to or has lasted longer than 8 weeks and the patient is in a suitable condition for the insertion of the tube.

The Degania™ Low Profile Gastrostomy Tube is a result of extensive expertise built through QMD's persistent search to merge clinical functionality with patients' needs. It is not just about innovation or technological features, it is about making life simpler, safer, and more comfortable for those who need to use enteral feeding products.

The Degania™ line of Low Profile Gastrostomy Tubes and accessories includes a comprehensive range of products that meet feeding, hydration, and venting requirements, allowing patients to continue daily life activities with minimum impact.

Our standard and special solutions for adults and children, the use of different materials, and the complete range of sizes, provide individual solutions for patients' needs.

### Degania™ G Button - Designed for support and performance



**1 Low Profile** Sits at skin level. Provides better comfort and mobility and a decreased risk of tube dislodgement.

**2 Apple-shaped balloon design** Helps reduce leakage by providing more surface contact with the stomach lining.

**3 Proximal anti reflux valve** Prevents the backflow of gastric fluid.

**4 Radiopaque striping** Allows visualization during initial placement.

**5 Sealing Cap** to close the access port.

**6 Recessed atraumatic tip** Prevents irritation of the stomach wall and allows the balloon to retract into the shaft of the tube, reducing discomfort during removal.

**Tapered distal tip** Enables easy, atraumatic insertion without a stiffener.



## INCREASED BALLOON DURABILITY

Tests performed on the Degania™ Low Profile Gastrostomy Tube are based on the ISO 20695-2020 standard. The tests showed considerably longer balloon durability in controlled laboratory conditions of simulated gastric fluid (pH 1.2 +/- 0.1).

### LPGT/B 12 Fr

First balloon bursts after 87 days


### LPGT/B 16 Fr

First balloon bursts after 83 days

### LPGT/B 24 Fr

First balloon bursts after 87 days

Surviving in simulated gastric solution for 83 days is equivalent to an actual lifetime of 332 days, and 87 days is equivalent to an actual lifetime of 348 days.



QMD is committed to enhancing patient safety by offering ISO 80369-3 compliant (ENFit®) products to reduce the risk of enteral administration misconnection also supported by the United States Food and Drug Administration (FDA), professional organizations, and clinicians around the world.



### Standardized Connectors

ISO 80369-3 standardized connectors specifically for enteral feeding.



### Reduce Misconnections

Special connector design ensures that enteral administration devices can connect only to enteral access devices.



### Reduce Accidental Disconnection

The locking threads help prevent the connector from accidentally detaching from the tubing.

## KIT CONCEPT

The Degania™ Percutaneous Replacement Gastrostomy Tube (G Button) Kit provides patients with access to everything needed for quick and easy insertion, with components that are compatible with one another.

### Degania™ Percutaneous Replacement Gastrostomy Tube (G Button) Kit includes:

- **Low Profile Gastrostomy Tube (G Button)**
- **Lubricating Gel Pack**
- **Gauze Pack**
- **Dual ENFit® Feeding Set with right angled secure lock connector**
- **Single ENFit® Feeding Set with straight secure lock connector**
- **ENFit® Syringe – 35cc**
- **Luer Slip Syringe – 10cc**



## ACCESSORIES

### Stoma Measuring Device



The Degania™ Stoma Measuring Device includes an external plate that allows you to determine the proper length of a well-established gastrostomy stoma tract to assist the proper selection of the Gastrostomy Button, thus minimizing leakage and associated complications.

### Extension Sets

Enteral feeding can be administered through Low Profile Gastrostomy Tubes using various types of Degania™ Extension Sets:

**INTERMITTENT GRAVITY DRIP** - the feed solution is administered over a set volume over 30-60 minutes several times a day, often by a given set only.

**CONTINUOUS** - the feed solution is administered over a period of 8-24 hours, often using an enteral feeding pump set to a prescribed rate.

**BOLUS** - the feed solution is administered over a 15-20-minute period, often via a syringe, several times a day; frequently in 150-200ml each session.

These methods can be used in combination and allow increased flexibility with the feeding regime to meet the requirements of the patient.



Degania™ Extension Feeding Sets are manufactured from high grade polyurethane tubing and are available with two types of secure-lock connectors and a tubing clamp. The sets are compatible with the most used Low Profile Gastrostomy tubes in the market.

# PRODUCT CATALOG

## Degania™ Low Profile Gastrostomy Tube - G Button - Non-Sterile

Product Number	Product Description	Size	Balloon Size	Length	Sterile	Inner bag	Carton Qty.
103302011263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	0.8cm	No	50 units	100 units
103302011463AQ		14FR	4ml		No	50 units	100 units
103302011663AQ		16FR	5ml		No	50 units	100 units
103302011863AQ		18FR	5ml		No	50 units	100 units
103302012063AQ		20FR	8ml		No	50 units	100 units
103302012463AQ		24FR	8ml		No	50 units	100 units
103302021263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	1.0cm	No	50 units	100 units
103302021463AQ		14FR	4ml		No	50 units	100 units
103302021663AQ		16FR	5ml		No	50 units	100 units
103302021863AQ		18FR	5ml		No	50 units	100 units
103302022063AQ		20FR	8ml		No	50 units	100 units
103302022463AQ		24FR	8ml		No	50 units	100 units
103302031263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	1.2cm	No	50 units	100 units
103302031463AQ		14FR	4ml		No	50 units	100 units
103302031663AQ		16FR	5ml		No	50 units	100 units
103302031863AQ		18FR	5ml		No	50 units	100 units
103302032063AQ		20FR	8ml		No	50 units	100 units
103302032463AQ		24FR	8ml		No	50 units	100 units
103302041263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	1.5cm	No	50 units	100 units
103302041463AQ		14FR	4ml		No	50 units	100 units
103302041663AQ		16FR	5ml		No	50 units	100 units
103302041863AQ		18FR	5ml		No	50 units	100 units
103302042063AQ		20FR	8ml		No	50 units	100 units
103302042463AQ		24FR	8ml		No	50 units	100 units
103302051263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	1.7cm	No	50 units	100 units
103302051463AQ		14FR	4ml		No	50 units	100 units
103302051663AQ		16FR	5ml		No	50 units	100 units
103302051863AQ		18FR	5ml		No	50 units	100 units
103302052063AQ		20FR	8ml		No	50 units	100 units
103302052463AQ		24FR	8ml		No	50 units	100 units
103302061263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	2.0cm	No	50 units	100 units
103302061463AQ		14FR	4ml		No	50 units	100 units
103302061663AQ		16FR	5ml		No	50 units	100 units
103302061863AQ		18FR	5ml		No	50 units	100 units
103302062063AQ		20FR	8ml		No	50 units	100 units
103302062463AQ		24FR	8ml		No	50 units	100 units

Degania™ Low Profile Gastrostomy Tube - G Button - Non-Sterile

Product Number	Product Description	Size	Balloon Size	Length	Sterile	Inner bag	Carton Qty.
103302071263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	2.3cm	No	50 units	100 units
103302071463AQ		14FR	4ml		No	50 units	100 units
103302071663AQ		16FR	5ml		No	50 units	100 units
103302071863AQ		18FR	5ml		No	50 units	100 units
103302072063AQ		20FR	8ml		No	50 units	100 units
103302072463AQ		24FR	8ml		No	50 units	100 units
103302081263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	2.5cm	No	50 units	100 units
103302081463AQ		14FR	4ml		No	50 units	100 units
103302081663AQ		16FR	5ml		No	50 units	100 units
103302081863AQ		18FR	5ml		No	50 units	100 units
103302082063AQ		20FR	8ml		No	50 units	100 units
103302082463AQ		24FR	8ml		No	50 units	100 units
103302091263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	2.7cm	No	50 units	100 units
103302091463AQ		14FR	4ml		No	50 units	100 units
103302091663AQ		16FR	5ml		No	50 units	100 units
103302091863AQ		18FR	5ml		No	50 units	100 units
103302092063AQ		20FR	8ml		No	50 units	100 units
103302092463AQ		24FR	8ml		No	50 units	100 units
103302101263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	3.0cm	No	50 units	100 units
103302101463AQ		14FR	4ml		No	50 units	100 units
103302101663AQ		16FR	5ml		No	50 units	100 units
103302101863AQ		18FR	5ml		No	50 units	100 units
103302102063AQ		20FR	8ml		No	50 units	100 units
103302102463AQ		24FR	8ml		No	50 units	100 units
103302111263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	3.5cm	No	50 units	100 units
103302111463AQ		14FR	4ml		No	50 units	100 units
103302111663AQ		16FR	5ml		No	50 units	100 units
103302111863AQ		18FR	5ml		No	50 units	100 units
103302112063AQ		20FR	8ml		No	50 units	100 units
103302112463AQ		24FR	8ml		No	50 units	100 units
103302121263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	4.0cm	No	50 units	100 units
103302121463AQ		14FR	4ml		No	50 units	100 units
103302121663AQ		16FR	5ml		No	50 units	100 units
103302121863AQ		18FR	5ml		No	50 units	100 units
103302122063AQ		20FR	8ml		No	50 units	100 units
103302122463AQ		24FR	8ml		No	50 units	100 units
103302131263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	4.5cm	No	50 units	100 units
103302131463AQ		14FR	4ml		No	50 units	100 units
103302131663AQ		16FR	5ml		No	50 units	100 units
103302131863AQ		18FR	5ml		No	50 units	100 units
103302132063AQ		20FR	8ml		No	50 units	100 units
103302132463AQ		24FR	8ml		No	50 units	100 units

**Degania™ Low Profile Gastrostomy Tube - G Button – Non-Sterile**

Product Number	Product Description	Size	Balloon Size	Length	Sterile	Inner bag	Carton Qty.
103302141263AQ	Low Profile Gastrostomy Tube - G Button	12FR	3ml	5.0cm	No	50 units	100 units
103302141463AQ		14FR	4ml		No	50 units	100 units
103302141663AQ		16FR	5ml		No	50 units	100 units
103302141863AQ		18FR	5ml		No	50 units	100 units
103302142063AQ		20FR	8ml		No	50 units	100 units
103302142463AQ		24FR	8ml		No	50 units	100 units
103302151463AQ	Low Profile Gastrostomy Tube - G Button	14FR	4ml	5.5cm	No	50 units	100 units
103302151863AQ		18FR	5ml		No	50 units	100 units

**Degania™ Percutaneous Replacement Gastrostomy Tube (G Button) Kit - Sterile**

Product Number	Product Description	Size	Length	ENFit®	Sterile	Box Qty.	Carton Qty.
103302011288AQ	Percutaneous Replacement Gastrostomy Tube (G Button) Kit	12FR	0.8cm	N/A	Yes	1 unit	10 units
103302021288AQ			1.0cm	N/A	Yes	1 unit	10 units
103302031288AQ			1.2cm	N/A	Yes	1 unit	10 units
103302041288AQ			1.5cm	N/A	Yes	1 unit	10 units
103302051288AQ			1.7cm	N/A	Yes	1 unit	10 units
103302061288AQ			2.0cm	N/A	Yes	1 unit	10 units
103302071288AQ			2.3cm	N/A	Yes	1 unit	10 units
103302081288AQ			2.5cm	N/A	Yes	1 unit	10 units
103302091288AQ			2.7cm	N/A	Yes	1 unit	10 units
103302101288AQ			3.0cm	N/A	Yes	1 unit	10 units
103302111288AQ			3.5cm	N/A	Yes	1 unit	10 units
103302121288AQ			4.0cm	N/A	Yes	1 unit	10 units
103302131288AQ			4.5cm	N/A	Yes	1 unit	10 units
103302141288AQ			5.0cm	N/A	Yes	1 unit	10 units
103302011488AQ	Percutaneous Replacement Gastrostomy Tube (G Button) Kit Balloon Standard Kit	14FR	0.8cm	N/A	Yes	1 unit	10 units
103302021488AQ			1.0cm	N/A	Yes	1 unit	10 units
103302031488AQ			1.2cm	N/A	Yes	1 unit	10 units
103302041488AQ			1.5cm	N/A	Yes	1 unit	10 units
103302051488AQ			1.7cm	N/A	Yes	1 unit	10 units
103302061488AQ			2.0cm	N/A	Yes	1 unit	10 units
103302071488AQ			2.3cm	N/A	Yes	1 unit	10 units
103302081488AQ			2.5cm	N/A	Yes	1 unit	10 units
103302091488AQ			2.7cm	N/A	Yes	1 unit	10 units
103302101488AQ			3.0cm	N/A	Yes	1 unit	10 units
103302111488AQ			3.5cm	N/A	Yes	1 unit	10 units
103302121488AQ			4.0cm	N/A	Yes	1 unit	10 units
103302131488AQ			4.5cm	N/A	Yes	1 unit	10 units
103302141488AQ			5.0cm	N/A	Yes	1 unit	10 units
103302151488AQ			5.5cm	N/A	Yes	1 unit	10 units
103302161488AQ			6.0cm	N/A	Yes	1 unit	10 units
103302171488AQ	6.5cm	N/A	Yes	1 unit	10 units		



Degania™ Percutaneous Replacement Gastrostomy Tube (G Button) Kit - Sterile

Product Number	Product Description	Size	Length	ENFit®	Sterile	Box Qty.	Carton Qty.
103302011888AQ	Percutaneous Replacement Gastrostomy Tube (G Button) Kit	18FR	0.8cm	N/A	Yes	1 unit	10 units
103302021888AQ			1.0cm	N/A	Yes	1 unit	10 units
103302031888AQ			1.2cm	N/A	Yes	1 unit	10 units
103302041888AQ			1.5cm	N/A	Yes	1 unit	10 units
103302051888AQ			1.7cm	N/A	Yes	1 unit	10 units
103302061888AQ			2.0cm	N/A	Yes	1 unit	10 units
103302071888AQ			2.3cm	N/A	Yes	1 unit	10 units
103302081888AQ			2.5cm	N/A	Yes	1 unit	10 units
103302091888AQ			2.7cm	N/A	Yes	1 unit	10 units
103302101888AQ			3.0cm	N/A	Yes	1 unit	10 units
103302111888AQ			3.5cm	N/A	Yes	1 unit	10 units
103302121888AQ			4.0cm	N/A	Yes	1 unit	10 units
103302131888AQ			4.5cm	N/A	Yes	1 unit	10 units
103302141888AQ			5.0cm	N/A	Yes	1 unit	10 units
103302151888AQ			5.5cm	N/A	Yes	1 unit	10 units
103302161888AQ			6.0cm	N/A	Yes	1 unit	10 units
103302012081AQ	Percutaneous Replacement Gastrostomy Tube (G Button) Kit	20FR	0.8cm	N/A	Yes	1 unit	10 units
103302022081AQ			1.0cm	N/A	Yes	1 unit	10 units
103302032081AQ			1.2cm	N/A	Yes	1 unit	10 units
103302042081AQ			1.5cm	N/A	Yes	1 unit	10 units
103302052081AQ			1.7cm	N/A	Yes	1 unit	10 units
103302062081AQ			2.0cm	N/A	Yes	1 unit	10 units
103302072088AQ			2.3cm	N/A	Yes	1 unit	10 units
103302082088AQ			2.5cm	N/A	Yes	1 unit	10 units
103302092081AQ			2.7cm	N/A	Yes	1 unit	10 units
103302102081AQ			3.0cm	N/A	Yes	1 unit	10 units
103302112081AQ			3.5cm	N/A	Yes	1 unit	10 units
103302122081AQ			4.0cm	N/A	Yes	1 unit	10 units
103302132081AQ			4.5cm	N/A	Yes	1 unit	10 units
103302142081AQ			5.0cm	N/A	Yes	1 unit	10 units
103302012481AQ	Percutaneous Replacement Gastrostomy Tube (G Button) Kit	24FR	0.8cm	N/A	Yes	1 unit	10 units
103302022481AQ			1.0cm	N/A	Yes	1 unit	10 units
103302032481AQ			1.2cm	N/A	Yes	1 unit	10 units
103302042481AQ			1.5cm	N/A	Yes	1 unit	10 units
103302052481AQ			1.7cm	N/A	Yes	1 unit	10 units
103302062481AQ			2.0cm	N/A	Yes	1 unit	10 units
103302072481AQ			2.3cm	N/A	Yes	1 unit	10 units
103302082481AQ			2.5cm	N/A	Yes	1 unit	10 units
103302092481AQ			2.7cm	N/A	Yes	1 unit	10 units
103302102481AQ			3.0cm	N/A	Yes	1 unit	10 units
103302124381AQ			3.5cm	N/A	Yes	1 unit	10 units
103302122481AQ			4.0cm	N/A	Yes	1 unit	10 units
103302132481AQ			4.5cm	N/A	Yes	1 unit	10 units
103302142481AQ			5.0cm	N/A	Yes	1 unit	10 units

DEGANIA™ Stoma Measuring Device - Sterile

Product Number	Product Description	Size	Balloon Size	Length	ENFit®	Sterile	Box Qty.	Carton Qty.
102148101280AQ	Stoma Measuring Device	12FR	5ml	205mm	No	Yes	2 units	100 units

DEGANIA™ Extension Feeding Sets - Sterile

Product Number	Product Description	Length	Secure Lock Connector	Clamp	Sterile	Box Qty.	Carton Qty.
1032430080AQ	Single ENFit® Port Extension Feeding Set	300mm	Straight	Yes	Yes	10 units	100 units
1032431080AQ			Right Angled	Yes	Yes	10 units	100 units
1032460080AQ		600mm	Straight	Yes	Yes	10 units	100 units
1032461080AQ			Right Angled	Yes	Yes	10 units	100 units
1032530080AQ	Dual (Y type) ENFit® Port Extension Feeding Set	300mm	Straight	Yes	Yes	10 units	100 units
1032531080AQ			Right Angled	Yes	Yes	10 units	100 units
1032560080AQ		600mm	Straight	Yes	Yes	10 units	100 units
1032561080AQ			Right Angled	Yes	Yes	10 units	100 units
1032030080AQ	Male Luer Lock Connector Port Extension Feeding Set	300mm	Straight	Yes	Yes	10 units	100 units
1032031080AQ			Right Angled	Yes	Yes	10 units	100 units
1032060080AQ		600mm	Straight	Yes	Yes	10 units	100 units
1032061080AQ			Right Angled	Yes	Yes	10 units	100 units
1032230080AQ			Straight	Yes	Yes	10 units	100 units
1032260080AQ			Straight	Yes	Yes	10 units	100 units
1032930080AQ	Bolus Universal Port Extension Feeding Set	300mm	Straight	Yes	Yes	10 units	100 units
1032960080AQ		600mm	Straight	Yes	Yes	10 units	100 units
1032730080AQ	2 Port (Y type) Continuous Extension Feeding Set	300mm	Straight	Yes	Yes	10 units	100 units
1032731080AQ			Right Angled	Yes	Yes	10 units	100 units
1032760080AQ		600mm	Straight	Yes	Yes	10 units	100 units
1032761080AQ			Right Angled	Yes	Yes	10 units	100 units
1032130080AQ	EN Lock Connector Port Extension Feeding Set	300mm	Straight	Yes	Yes	10 units	100 units
1032131080AQ			Right Angled	Yes	Yes	10 units	100 units
1032160080AQ		600mm	Straight	Yes	Yes	10 units	100 units
1032161080AQ			Right Angled	Yes	Yes	10 units	100 units

