


KIMAL
delivering healthcare innovation



YOUR ALTIUS CVC RANGE



REDUCE
PATIENT **RISKS**
WITH ALTIUS



WHY TO CHOOSE ALTIUS™



We have taken the feedback from health care professionals and designed **Altius** - the **first innovation in CVCs in over 40 years**.

Altius is manufactured utilising our proprietary multitube technique to achieve perfectly round lumens and therefore achieving laminar flow.

This allows for higher flows than conventional extruded catheters especially in 4, 5 & 7 lumen configuration.

Altius Proactiv+ is a chlorhexidine free non eluting covalently bonded antimicrobial to **protect your patient for up to 30 days**.

Altius Proactiv+™ unique covalent bonding ensures no systemic release into patient's bloodstream and ensures sustained antimicrobial activity.

THE LUERSAFE CONNECTORS


Integrated connector to prevent the risk of air embolism and infection during manipulation of the lines. No need to clamp and unclamp the lines.

COLOUR CODED LUMENS

To easily identify what drug goes through what lumen and reduce to a minimum the manipulation of the lines.

MULTITUBE TECHNOLOGY

Perfectly round lumens to deliver higher flows - over 21% higher flows compared to leading competitor in 5 lumen*.



Altius Proactiv+ combines a non thrombogenic (PEG) with the antimicrobial polyhexamethyl biguanide (PHMB) to provide superior antimicrobial effectiveness.

WHY PHMB

It is a fast acting, broad spectrum biocide effective at low concentrations against gram positive and gram negative bacteria, fungi and viruses and widely used in antimicrobial medical devices.

HOW DOES PROACTIV+ WORK?


PHMB is co-polymerised with a non thrombogenic PEG (polyethylene glycol) to prevent platelet adhesion reducing the opportunity for bacterial growth. It reduces biofilm formation and minimises fibrin sheath.

The patented non-eluting crosslinked covalent bonding ensures patient protection for the life of the catheter.

THE RANGE

- **Adults:** 1 to 5 and 7 lumen.
- **Paediatric:** 1 to 5 lumen.

For french size, length and flow rates please check ordering guide.



THE RANGE

- **Adults:** 1 to 5 and 7 lumen.

For french size, length and flow rates please check ordering guide.



SUPERIOR INSERTION KIT

The superior insertion kit has been designed to reduce the risks during insertion.




Super Sharp Echogenic Needle: Easy visualisation under the ultrasound to prevent posterior wall puncture.



5m Guidewire Introducer: Luerslip syringe for closed system.



Altius Guidewire Dispenser: Ergonomic design to facilitate guidewire advancement even with wet gloves.



Kink Resistant Nitinol Guidewire: To prevent guidewire kinking even in difficult patient's anatomy.



N11 Mini Scalpel: To nick skin prior to dilation.




Hydrophilic Dilator: For smooth dilation to prevent insertion trauma.



Secondary Fixation Wing

*comparison of market leader 5 lumen 8.5 french 16cm published data with Altius Classic 5 lumen 8.5Fr.

YOUR ALTIUS HP CVC RANGE



REDUCE
PATIENT **RISKS**
WITH ALTIUS



WHY TO CHOOSE ALTIUS™ HIGH PRESSURE



Altius HP has all the benefits of Altius with a **High Pressure purple line** to deliver contrast media and fluids at high pressures of upto **300 psi** with a **flow of 10ml per second**.

Altius HP can deliver contrast media up to a viscosity rating of 11.8cP units.

Altius HP is the first CVC with a dedicated purple line compatible to deliver contrast media, reducing the need for additional cannulations in critical patients that are hypovolemic and difficult to cannulate.


Altius HP purple line makes it very easy to identify across functional teams which lumen is the contrast media compatible lumen.

Altius HP is manufactured utilising our proprietary multitube technique to achieve perfectly round lumens and therefore laminar flow. This allows for higher flows especially in 4 & 5 lumen configurations.

All the benefits of **Altius Proactiv+** and a **High Pressure purple line** to deliver contrast media and fluids at high pressures.

Proactiv+ is a chlorhexidine free, non-eluting, covalently bonded antimicrobial to protect your patient for up to 30 days.

Proactiv+™ unique covalent bonding ensures no systemic release into patient's bloodstream.



HIGH PRESSURE LINE

A high pressure line to deliver power injectable contrast media up to 300psi and a maximum of 10ml/sec.

THE LUERSAFE CONNECTORS

Integrated connector to prevent the risk of air embolism and infection during manipulation of the lines. No need to clamp and unclamp the lines.

COLOUR CODED LUMENS

To easily identify what drug goes through what lumen and reduce to a minimum the manipulation of the lines.

MULTITUBE TECHNOLOGY

Perfectly round lumens to deliver higher flows - over 21% higher flows compared to leading competitor in 5 lumen*.

Proactiv+ combines a non thrombogenic (PEG) with the antimicrobial polyhexamethyl biguanide (PHMB) to provide superior antimicrobial effectiveness.

WHY PHMB

It is a fast acting, broad spectrum biocide effective at low concentrations against gram positive and gram negative bacteria, fungi and viruses and widely used in antimicrobial medical devices.

HOW DOES PROACTIV+ WORK?

PHMB is co-polymerised with a non thrombogenic PEG (polyethylene glycol) to prevent platelet adhesion reducing the opportunity for bacterial growth. It reduces biofilm formation and minimises fibrin sheath.

The patented non-eluting crosslinked covalent bonding ensures patient protection for the life of the catheter.

THE RANGE

• **Adults:** 1 to 5 lumen.

For french size, length and flow rates please check ordering guide.



*comparison of market leader 5 lumen 8.5 french 16cm published data with Altius Classic 5 lumen 8.5Fr.

THE RANGE

• **Adults:** 1 to 5 lumen.

For french size, length and flow rates please check ordering guide.



SUPERIOR INSERTION KIT

The superior insertion kit has been designed to reduce the risks during insertion.



Super Sharp Echogenic Needle: Easy visualisation under the ultrasound to prevent posterior wall puncture.



5m Guidewire Introducer: Luerslip syringe for closed system.



Altius Guidewire Dispenser: Ergonomic design to facilitate guidewire advancement even with wet gloves.



Kink Resistant Nitinol Guidewire: To prevent guidewire kinking even in difficult patient's anatomy.



N11 Mini Scalpel: To nick skin prior to dilation.



Hydrophilic Dilator: For smooth dilation to prevent insertion trauma.



Secondary Fixation Wing